

Cardiac Science Powerheart® G5 Automated External Defibrillator (AED)

What happens in the first few minutes following a sudden cardiac arrest (SCA) can make the difference between life and death. The Powerheart® G5 is the first AED that combines variable escalating energy with fast shock times and real-time CPR feedback. Easy, intuitive operation and user-paced prompts enable responders in workplaces and public spaces to act quickly and confidently to give an SCA victim the best chance of a favorable outcome.

Features and benefits

Easy to use

Users can successfully perform each critical step of rescue with the Powerheart G5's real-time CPR feedback, user-paced instruction, intuitive pad design and automatic shock delivery. (Semi-automatic model also available).

Rugged and reliable

Ready to go whenever you need it with a durable design tested to military standards and an IP55 rating for protection from dust and water. The highly visible Rescue Ready® indicator confirms your G5 has passed rigorous daily, weekly, and monthly self-tests.

Customisable

Customise the prompts, shock protocols, and CPR feedback to match user skills and preferences and easily update these settings if Resuscitation Guidelines change.

Incredibly smart

Rely on the Powerheart G5's advanced technology with customised escalating energy and fast shock times. Easily switch to an alternate language* during a rescue with the push of a button.

Warranty

The G5 is built to last with a seven-year device warranty and fouryear full replacement warranty on the medical grade battery.

*Available in selected languages



To find out more contact **0161 930 6010** or visit **www.boclifeline.co.uk**

Powerheart G5 AED Technical specifications

Defibrillator	
Operations	Fully automatic and semi-automatic versions
Waveform	STAR® biphasic truncated exponential
Therapy range	95J to 354J (adult)
Factory default (nominal)	Escalating Variable Energy (VE): 200VE, 200VE
Shock delivery	Within 10 seconds (typically) post-CPR
Voice prompts	RescueCoach™ prompts provide user-paced instruction
Text screen	Displays rescue prompts and critical rescue information
Audible alerts	Voice prompt, system alert
Synchronized shock	Built-in automatic synchronisation feature
Pacemaker pulse detection	Yes
Customisable	Voice prompt, CPR settings, and shock protocols via AED Manager
CPR feedback	Real-time guidance for CPR rate/ depth (optional), metronome for CPR rate
Visible indicators	Rescue Ready®, battery, service, and pad status
Warranty	7 years

Physical/environmental

Dimensions (L x W x D)	30 cm x 23 cm x 9 cm
	(11.8 in x 9.1 in x 3.5 in)
Weight	2.5 kg (5.5 lbs), including battery
	and pads
Dust/Water resistance	IP55 (IEC 60529)
Operating temperature	0°C to +50°C (32°F to +122°F)
Altitude	-382 m (-1,253 ft) to
	4,594 m (15,073 ft)
Drop	MIL-STD-810G, 516.6 (Procedure IV)
	for 1.22 m drop
Shocks	MIL-STD-810G, 516.6 (Procedure 1)
Vibration (sine and random)	MIL-STD-810G, 514.6
	(Procedure 1, Cat 24)
EMI (radiated/immunity)	IEC 60601-2-4

P	a	d	S

· ·	Disposable, non-polarized (pads can
	be placed in either position)
Shelf life	2 years

Battery

Туре	Intellisense [®] lithium battery
Guarantee	4-year, full operational replacement
Capacity	500 shocks at 150VE; 450 shocks
	at 200VE

Automatic self-tests

Daily	Battery, pads (presence and function), electrical circuitry, and software
Weekly	Partial energy charge
Monthly	Full energy charge cycle

Event documentation

Internal memory	90 minutes of rescue data, multiple
	rescue functionality
ECG and rescue review	Viewable via AED Manager software
Communications	USB cable or USB drive

Each AED package includes (1) defibrillator, (1) Intellisense battery, (1) pair of adult defibrillation pads**, (1) AED Manager software, (1) AED manual, Getting Started card, and (1) Steps to a rescue card.

 $\begin{tabular}{c} \textbf{C} & \textbf{Cardiac Science, Powerheart, RescueCoach, Rescue Ready, Intellisense, and } \\ & \circ \circ \circ \circ & \textbf{STAR are trademarks of Cardiac Science Corporation. All Rights Reserved.} \\ \end{tabular}$

BOC Healthcare is an authorised distributor of Automated External Defibrillators manufactured by Cardiac Science.

BOC Healthcare, Great Britain & N. Ireland Customer Service Centre, Priestley Road, Worsley, Manchester M28 2UT, United Kingdom Tel 0161 930 6010 bochealthcare-uk@boc.com, www.boclifeline.co.uk BOC Healthcare, Republic of Ireland PO Box 201, Bluebell, Dublin 12, Republic of Ireland Tel 01 409 1827, Fax 01 409 1854 healthcareinfo.ie@boc.com, www.BOConline.ie

 $[\]ensuremath{^{\star\star}}\xspace$ Optional CPR device available at point of order.